  Microlife UAB
P. Lukšio g. 32
08222 Vilnius / Lithuania

 Microlife AG
Espenstrasse 139
9443 Widnau / Switzerland
www.microlife.com

CE0044

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- ① Measuring sensor / Removable cap
- ② START button
- ③ Display
- ④ ON/OFF button
- ⑤ Battery compartment cover
- ⑥ All segments displayed
- ⑦ Memory
- ⑧ Ready for measuring
- ⑨ Measurement complete
- ⑩ Surface Temperature Indication
- ⑪ Low battery indicator
- ⑫ Changing between Celsius and Fahrenheit
- ⑬ Recall mode
- ⑭ Recall the last 12 Readings
- ⑮ Measured temperature too high
- ⑯ Measured temperature too low
- ⑰ Ambient temperature too high
- ⑱ Ambient temperature too low
- ⑲ Error function display
- ⑳ Blank display
- ㉑ Flat battery
- ㉒ Replacing the battery

Read the instructions carefully before using this device.

Type BF applied part

Keep dry

This Microlife thermometer is a high quality product incorporating the latest technology and tested in accordance with international standards. With its unique technology, this device can provide a stable, heat-interference-free reading with each measurement. The device performs a self-test every time it is switched on to always guarantee the specified accuracy of any measurement.

This Microlife Thermometer is intended for the periodic measurement and monitoring of human body temperature. It is intended for use on people of all ages.

This thermometer has been clinically tested and proven to be safe and accurate when used in accordance to the operating instruction manual.

Please read through these instructions carefully in order for you to understand all functions and safety information.

Table of Contents

1. The Advantages of this Thermometer
2. Important Safety Instructions
3. How this Thermometer measures Temperature
4. Control Displays and Symbols
5. Changing between Forehead and Ear Mode
6. Directions for Use
7. Changing between Celsius and Fahrenheit
8. How to recall 12 readings in Memory Mode
9. Error Messages
10. Cleaning and Disinfecting
11. Battery Replacement
12. Guarantee
13. Technical Specifications
14. www.microlife.com
Guarantee Card (see Back Cover)

1. The Advantages of this Thermometer

Multiple uses (wide range of measurement)

This thermometer offers a wide measurement range feature from 0 - 100.0 °C / 32.0 - 212.0 °F, meaning the unit can be used as a forehead or an ear thermometer to measure body temperature or it can be used to measure surface temperature of the following:

- Milk surface temperature in a baby's bottle
- Surface temperature of a baby's bath
- Ambient temperature

Probe cover free

This thermometer is more user-friendly and more cost-effective since a probe cover is not required.

Measures in a matter of seconds

The innovative infrared technology allows measurement of ear temperature in only 1 second and forehead temperature in only 3 seconds.

Accurate and reliable

The unique probe assembly construction incorporates an advanced infrared sensor, ensuring that each measurement is accurate and reliable.

Gentle and easy to use

- The ergonomic design enables simple and easy use of the thermometer.
- This thermometer can even be used on a sleeping child without causing any interruption.
- This thermometer is quick, therefore child-friendly.

Auto-Display Memory

The last reading is automatically displayed for 2 seconds when the unit is switched ON.

Multiple readings recall

Users will be able to recall the last 12 readings when entering the recall mode, enabling efficient tracking of temperature variations.

Safe and hygienic

- No risk of broken glass or mercury ingestion.
- Completely safe for use on children.
- Cleaning of the probe can be done with an alcohol-moistened cotton tissue, making this thermometer completely hygienic for use by the whole family.

Fever alarm

10 short beeps and a red LCD backlight alert the patient that he/she may have a temperature equal to or higher than 37.5 °C.

2. Important Safety Instructions

- Follow instructions for use. This document provides important product operation and safety information regarding this device. Please read this document thoroughly before using the device and keep for future reference.
- This device may only be used for the purposes described in these instructions. The manufacturer cannot be held liable for damage caused by incorrect application.

- **Never immerse this device in water or other liquids. For cleaning please follow the instructions in the «Cleaning and Disinfecting» section.**
- Do not use this device if you think it is damaged or notice anything unusual.
- Never open this device.
- Earwax in the ear canal may cause a lower temperature reading. Therefore it is important to ensure the subject's ear canal is clean.
- A basic physiologic effect called vasoconstriction can occur in the early stages of fever, resulting in a cool skin effect, the recorded temperature using a forehead thermometer can, therefore, be unusually low.
- If the measurement result is not consistent with the patient's finding or unusually low, repeat the measurement every 15 minutes or double check the result by another core body temperature measurement.
- This device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the «Technical Specifications» section.
- Ensure that children do not use this device unsupervised; some parts are small enough to be swallowed. Be

aware of the risk of strangulation in case this device is supplied with cables or tubes.

- Do not use this device close to strong electromagnetic fields such as mobile telephones or radio installations. Keep a minimum distance of 3.3 m from such devices when using this device.
- Protect it from:
 - extreme temperatures
 - impact and dropping
 - contamination and dust
 - direct sunlight
 - heat and cold
- If the device is not going to be used for a prolonged period the battery should be removed.

WARNING: The measurement results given by this device is not a diagnosis. It is not replacing the need for the consultation of a physician, especially if not matching the patient's symptoms. Do not rely on the measurement result only, always consider other potentially occurring symptoms and the patient's feedback. Calling a doctor or an ambulance is advised if needed.

3. How this Thermometer measures Temperature

This thermometer measures forehead and ear temperature.

Forehead measurement

This thermometer measures infrared energy radiated from the forehead. This energy is collected through the lens and converted to a temperature value.

Ear measurement

This thermometer measures infrared energy radiated from the eardrum and the surrounding tissue. This energy is collected through the lens and converted to a temperature value. The measured reading obtained directly from the eardrum (Tympanic Membrane) ensures the most accurate ear temperature.

Measurements taken from the surrounding tissue of the ear canal generate lower readings and may result in misdiagnosis of fever.

To avoid an inaccurate measurement

1. Switch on the thermometer by pressing the ON/OFF button ④.
2. After one beep is heard (and the temperature scale icon is flashing), straighten the ear canal by gently pulling the middle of the ear back and up.
3. Place the probe ① firmly into the ear canal, press the START button ② and keep the probe in the ear until the thermometer beeps to signal the completion of the measurement.

4. Control Displays and Symbols

- **All segments displayed** ⑥: Press the ON/OFF button ④ to turn on the unit; all segments will be shown for 2 seconds.
- **Memory** ⑦: The last reading will automatically be shown on the display for 2 seconds.
- **Ready for measurement** ⑧: The unit is ready for measurement, the «°C» or «°F» icon will keep flashing while the mode icon (ear or forehead) will be displayed.
- **Measurement complete** ⑨: The reading will be shown on the display ③ with the «°C» or «°F» icon and the mode icon steady. The unit is ready for the next measurement as soon as the «°C» or «°F» icon is flashing again.
- **Surface temperature indication** ⑩: A crossed-ear icon will appear on the display ③ when the reading in the ear mode is outside the range 32.0 ~ 42.2 °C / 89.6 ~ 108.0 °F.
- **Low battery indicator** ⑪: When the unit is turned on, the battery icon will keep flashing to remind the user to replace the battery.

5. Changing between Forehead and Ear Mode

To change from forehead to ear mode, take off the removable cap ① from the thermometer tip. The device will immediately switch to ear mode. The ear icon will be shown

on the display and the «°C» or «°F» icon will flash to indicate that the device is ready for measurement.

To change back to the forehead mode, place the removable cap ① on the thermometer tip and the device will immediately switch back to forehead mode. The forehead icon will be shown on the display and the «°C» or «°F» icon will flash to indicate that the device is ready for measurement.

6. Directions for Use

1. Press the ON/OFF button ④. The display ③ is activated to show all segments for 2 seconds.
2. The last measurement reading will automatically be shown on the display for 2 seconds with the «M» icon ⑦.
3. When the «°C» or «°F» icon is flashing, a beep sound is heard and the thermometer is ready for measuring ⑧.

Forehead measurement

- ▶ Place the measuring sensor gently at the center of the forehead. Please remove any hair, sweat or dirt from the forehead before measuring to improve the accuracy of the readings.
- ▶ Press the START button ② once and release it. A long beep sound will emit after 3 seconds to acknowledge the completion of measurement.
- ▶ Read the recorded temperature from the LCD display.

NOTE:

- **Patients and thermometer should stay in similar room condition for at least 30 minutes.**
- Don't take a measurement while or immediately after nursing a baby.
- Don't use the thermometer in high humidity environments.
- Patients should not drink, eat or exercise before/while taking the measurement.
- Before placing the thermometer probe onto the measurement area, remove dirt, hair, or sweat.
- Don't move the measurement device from the measuring area before hearing the termination beep.
- Always take the temperature in the same location, since temperature readings may vary according to locations.
- 10 short beeps and a red LCD backlight alert the patient that he/she may have a temperature equal to or higher than 37.5 °C.

Ear measurement

- ▶ Straighten the ear canal by pulling the ear up and back to give a clear view of the eardrum.
 - For children under 1 year; pull the ear straight back.
 - For children of 1 year to adult; pull the ear up and back.
- Also refer to the short instructions at the front.

- ▶ While gently pulling the ear, insert the probe snugly into the ear canal (**about 1 sec.**) and **immediately** press the START button ②. Release the button and wait for the beep sound. This indicates the end of the measurement.
- ▶ Remove the thermometer from the ear canal. The display shows the measured temperature ⑨.

NOTE:

- **To ensure accurate readings, wait at least 30 sec. after 3-5 continuous measurements.**
- Accumulation of earwax on the probe can result in less accurate temperature readings or cross infection between users. **Therefore, it is essential that the probe is clean before each measurement.** For cleaning, follow the instructions in the «Cleaning and Disinfecting» section.
- **After cleaning the measuring sensor ① with alcohol, wait 5 minutes before taking the next measurement,** in order to allow the thermometer to reach its operating reference temperature.
- For an infant, it is best to have the child lying flat with his/her head sideways so the ear is facing upwards. For an older child or adult, it is best to stand behind and slightly to the side of the patient.
- Always take the temperature in the same ear, since the temperature readings may be different from ear to ear.

- In the following situations it is recommended that three temperatures in the same ear be taken and the highest one taken as the reading:
 1. New born infants in the first 100 days.
 2. Children under three years of age with a compromised immune system and for whom the presence or absence of fever is critical.
 3. When the user is learning how to use the thermometer for the first time until he/she has familiarized himself/herself with the device and obtains consistent readings.
 4. If the measurement is surprisingly low.

7. Changing between Celsius and Fahrenheit

This thermometer can display temperature readings in either Fahrenheit or Celsius. To switch the display between °C and °F, simply turn OFF the unit, **press and hold** the START button ② for 5 seconds; after 5 seconds, the current measurement scale («°C» or «°F» icon) will flash on the display ⑫. Change the measurement scale between °C and °F by pressing the START button ②. When the measurement scale has been chosen, wait for 5 seconds and the unit will automatically enter the «ready for measuring» mode.

8. How to recall 12 readings in Memory Mode

This thermometer can recall the last 12 readings.

- **Recall mode ⑬:** Press the START button ② to enter Recall mode when power is off. The memory icon «M» will flash.
- **Reading 1 - the last reading ⑭:** Press and release the START button ② to recall the last reading. Display 1 alone with memory icon.
- **Reading 12 - readings in succession:** Press and release the START button ② consecutively to recall the last 12 readings in succession.

Pressing and releasing the START button ② after the last 12 readings have been recalled will resume the above sequence from reading 1.

9. Error Messages

- **Measured temperature too high ⑮:** Displays «H» when measured temperature is higher than 100.0 °C / 212.0 °F in the ear mode or 42.2 °C / 108.0 °F in the forehead mode.
- **Measured temperature too low ⑯:** Displays «L» when measured temperature is lower than 0 °C / 32.0 °F in the ear mode or 34.0 °C / 93.2 °F in the forehead mode.

- **Ambient temperature too high** (17): Displays «H» in conjunction with the «▲» when ambient temperature is higher than 40.0 °C / 104.0 °F.
- **Ambient temperature too low** (18): Display «L» in conjunction with the «▼» when ambient temperature is lower than 10.0 °C / 50.0 °F in the ear mode or 16 °C / 60.8 °F in the forehead mode.
- **Error function display** (19): The system has a malfunction.
- **Blank display** (20): Check if the battery has been loaded correctly. Also check polarity (<+> and <->) of the battery.
- **Flat battery indication** (21): If the steady battery icon is the only symbol shown on the display, the battery should be replaced immediately.

10. Cleaning and Disinfecting

Use an alcohol swab or cotton tissue moistened with alcohol (70% Isopropyl) to clean the thermometer casing and the measuring sensor. Ensure that no liquid enters the interior of the device. Never use abrasive cleaning agents, thinners or benzene for cleaning and never immerse the device in water or other cleaning liquids. Take care not to scratch the surface of the sensor lens and the display.

11. Battery Replacement

This device is supplied with one lithium battery, type CR2032. Replace with a new CR2032 battery when the flashing battery symbol appears on the display (21). Remove the battery cover (22) by sliding it in the direction shown. Insert the new battery with the + at the top.

Batteries and electronic devices must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

12. Guarantee

This device is covered by a **2 year guarantee** from the date of purchase. During this guarantee period, at our discretion, Microlife will repair or replace the defective product free of charge.

Opening or altering the device invalidates the guarantee. The following items are excluded from the guarantee:

- Transport costs and risks of transport.
- Damage caused by incorrect application or non-compliance with the instructions for use.
- Damage caused by leaking batteries.
- Damage caused by accident or misuse.
- Packaging/storage material and instructions for use.

- Regular checks and maintenance (calibration).
- Accessories and wearing parts: Batteries, Probe cover (optional).

Should guarantee service be required, please contact the dealer from where the product was purchased, or your local Microlife service. You may contact your local Microlife service through our website:

www.microlife.com/support

Compensation is limited to the value of the product. The guarantee will be granted if the complete product is returned with the original invoice. Repair or replacement within guarantee does not prolong or renew the guarantee period. The legal claims and rights of consumers are not limited by this guarantee.

13. Technical Specifications

Type:	Dual Mode Thermometer IFR 100
Measurement range:	Ear Mode: 0-100.0 °C / 32-212.0 °F Forehead Mode: 34-42.2 °C / 93.2-108.0 °F
Resolution:	0.1 °C / °F
Measurement accuracy:	Laboratory: ±0.2 °C, 32.0 ~ 42.2 °C ±0.4 °F, 89.6 ~ 108.0 °F

Display: Liquid Crystal Display, 4 digits plus special icons

Acoustic: The unit is turned ON and ready for the measurement: 1 short beep.
Complete the measurement: 1 long beep (1 sec.) if the reading is less than 37.5 °C / 99.5 °F, 10 short «beep» sounds, if the reading is equal to or greater than 37.5 °C / 99.5 °F.
System error or malfunction: 3 short «bi» sounds.

Memory: Auto-Display the last measured temperature
12 readings recall in the Memory Mode

Backlight: The display light will be GREEN for 4 seconds, when the unit is turned ON.
The display light will be GREEN for 5 seconds, when a measurement is completed with a reading less than 37.5 °C / 99.5 °F.
The display light will be RED for 5 seconds, when a measurement is completed with a reading equal to or higher than 37.5 °C / 99.5 °F.

Operating conditions: Ear Mode: 10-40.0 °C / 50-104.0 °F
Forehead Mode: 16-40.0 °C / 60.8-104.0 °F

Storage conditions: -25 - +55 °C / -13 - +131 °F
15 - 95 % relative maximum humidity

Automatic Switch-off: Approx. 1 minute after last measurement has been taken.

Battery: 1 x CR2032 battery 3V

Battery lifetime: approx. 2000 measurements (using a new battery)

Dimensions: 107 x 50 x 34 mm

Weight: 55 g (with battery), 52 g (w/o battery)

IP Class: IP21

Reference to standards: EN 12470-5; ASTM E1965; IEC 60601-1; IEC 60601-1-2 (EMC); IEC 60601-1-11

Expected service life: 5 years or 12000 measurements

This device complies with the requirements of the Medical Device Directive 93/42/EEC.
Technical alterations reserved.
According to the Medical Product User Act a biennial technical inspection is recommended for professional users.
Please observe the applicable disposal regulations.

14. www.microlife.com

Detailed user information about our thermometers and blood pressure monitors as well as services can be found at www.microlife.com.

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